Baleen™ polyp trap
Piège à polypes Baleen™
Baleen™ Polypenfalle
Dispositivo per l'estrazione dei polipi Baleen™
Trampa de succión de pólipos Baleen™
Captador de pólipos Baleen™
Baleen™-polyropsamling
Baleen™-poliepvangen
Baleen™ polip kapanı
Baleen™ 응중 트랩
Экстрактор полипов Baleen™
Baleen™ ポリープトラップ

Reorder No. NBPT520
Référence de commande NBPT520
Nachbestell-Nr. NBPT520
N. di riordino NBPT520
N.º de pedido NBPT520
N.º de encomenda NBPT520
Genbestillingsnr. NBPT520
Nabestelnr. NBPT520
Yeni Sipariş No. NBPT520
주문 번호 NBPT520
Номер для повторного заказа NBPT520
再注文番号 NBPT520

INSTRUCTIONS FOR USE

MODE D'EMPLOI

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BRUGSANVISNING

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KULLANIM TALIMATLARI

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email global@usendoscopy.com

www.usendoscopy.com

This device is not made with natural rubber latex.

Intended use:

The single use Baleen™ polyp trap is used for suction retrieval of endoscopically removed polyps.

Warnings and Precautions:

- Consult the medical literature relative to techniques, complications and hazards prior to the performance of any endoscopic procedure.
- Endoscopic procedures should only be performed by persons having adequate training and familiarity with endoscopic techniques and polyp retrieval.
- Polyps less than or equal to 2mm may get suctioned through the mesh.
- Do not attempt to reuse, reprocess, refurbish, remanufacture or resterilize this device. US Endoscopy did not design this device nor is it intended to be reused, reprocessed, refurbished, remanufactured, or resterilized. Performing such activities on this disposable medical device presents a safety risk to patients (i.e. compromised device integrity, cross-contamination, infection).
- Inspect the package for shipping or handling damage. If damage is evident, do not use this device, save it for return, and contact your local Product Specialist.
- When removing a specimen from the trap, make sure that the entire net is grasped prior to removing trap from suction tubing to prevent the loss of specimen.

Contraindications:

- Contraindications include, but are not limited to, those specific to any endoscopic procedure.

Prior to Use:

1. Inspect and familiarize yourself with the device.
2. Read the Instructions for Use.
3. If there is evidence of damage (i.e. deformed, ripped, missing component, or damaged packaging), **do not use this product and contact your local Product Specialist.**

Directions for Use:

1. Install the Baleen™ polyp trap just prior to suctioning a polyp through the endoscope.
2. Remove the Baleen™ polyp trap from the holding tray.
3. Disconnect the suction tubing from the suction connector (barb) located on the umbilical section of the endoscope.
4. Insert the narrow tip end of the Baleen™ polyp trap into the suction tubing until you meet resistance (See Fig. 1).
5. Fold the excess netting protruding from the end of the suction tubing back over the end of the suction tubing (See Fig 2).
6. Grab the flange-end of the Baleen™ polyp trap applicator and remove the applicator from the end of the suction tubing leaving the netting in place (See Fig 3).
7. Discard the applicator after use.
8. Attach the suction tubing with the intact Baleen™ polyp trap to the suction connector (barb) of the endoscope (See Fig 4).
9. The Baleen™ polyp trap is now in place and ready for use.
10. After the polyp has been suctioned into the Baleen™ trap, disconnect the suction tubing from the umbilicus of the endoscope (holding onto the netting) (See Fig. 5).
11. To remove the polyp, grasp **all** the exposed netting and pull the netting out of the suction tubing.
Caution: Not grasping all the netting could potentially allow the polyp to fall out of netting.
12. Reattach suction line to the suction connector (barb) on the endoscope.
13. Check to confirm polyp is retrieved by carefully opening the netting.
14. Place the polyp (and netting-if allowed per facility guidelines) into the pathology specimen container.
15. Specimens should be prepared for histologic evaluation per institutional guidelines.



Figure 1


Figure 2

Figure 3


Figure 4

Figure 5

Product Disposal:

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Issued Date: February 2018

Warning:

An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact US Endoscopy to determine if additional information is available.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

US Endoscopy, a wholly owned subsidiary of STERIS Corporation.

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Listening...and delivering solutions®, le logo US Endoscopy listening...and delivering solutions, ainsi que toutes les marques portant le signe ® ou ™ sont des marques déposées auprès de l'U.S. Patent and Trademark Office ou sont des marques commerciales détenues par US Endoscopy.
Listening...and delivering solutions®, das US Endoscopy listening...and delivering solutions-Design und alle mit ® oder ™ bezeichneten Marken sind beim U.S. Patent und Trademark Office (US-Patent- und Markenamt) eingetragene Marken oder Warenzeichen, die sich im Besitz von US Endoscopy befinden.
Listening...and delivering solutions®, o design US Endoscopy listening...and delivering solutions e todas as marcas assinaladas com ® ou ™ estão registadas no Departamento de Registos de Patentes dos EUA (U.S. Patent and Trademark Office) ou são marcas comerciais detidas pela US Endoscopy.
Listening...and delivering solutions®, US Endoscopy listening...and delivering solutions-designet og alle varemærker markeret med ® eller ™ er registreret i USA's patentregister (U.S. Patent and Trademark Office) eller er varemærker, der tilhører US Endoscopy.
Listening...and delivering solutions®, het ontwerp Listening...and delivering solutions van US Endoscopy en alle aanduidingen met ® of ™ zijn gedeponeerd bij het Patent and Trademark Office in de Verenigde Staten of zijn handelsmerken van US Endoscopy.
Listening...and delivering solutions®, US Endoscopy listening...and delivering solutions 디자인, 그리고® 또는 ™으로 표시된 모든 상표는 미국 특허청 (U.S. Patent and Trademark Office)에 등록되어 있거나 US Endoscopy가 소유하고 있는 상표입니다.