Biopsy Forceps

Biopsy Forceps

Pince à biopsie **Reorder No. 00711201** (230cm)
 Biopsiezange **Reorder No. 00711202** (230cm)
 Pinza per biopsia **Reorder No. 00711203** (160cm)
 Fórceps para biopsia **Reorder No. 00711204** (160cm)
 Pinça de biopsia
 Biopsitang
 Biopsietang
 Biyopsi Forsepsi
 생검 겸자
 Щипцы биопсийные одноразовые

Biopsy Forceps – enteroscope

Pince à biopsie - entéroscope **Reorder No. 00711291** (350cm)
 Biopsiezange - Enteroskop
 Pinza per biopsia - enteroscopio
 Fórceps para biopsia – enteroscopia
 Pinça de biopsia - enteroscópio
 Biopsitang - enteroskop
 Biopsietang - enteroscoop
 Biyopsi Forsepsi - enteroskop
 생검 겸자 - 내시경
 Щипцы биопсийные одноразовые для использования с энтероскопом

Biopsy Forceps – sheathed

Pince à biopsie - gainée **Reorder No. 00711205** (230cm)
 Biopsiezange – ummantelt
 Pinza per biopsia - rivestita
 Fórceps para biopsia - con protección
 Pinça de biopsia - revestida
 Biopsitang – sammenfoldet
 Biopsietang – bekleed
 Biyopsi Forsepsi – kılıflı
 생검 겸자 - 초상
 Щипцы биопсийные в защитной трубке

Biopsy Forceps – hot

Biopsy Forceps – hot

Pince à biopsie – chaude **Reorder No. 00711211** (230cm)
 Biopsiezange - für Hochfrequenzstrom
 Pinza per biopsia - a caldo
 Fórceps para biopsia – caliente
 Pinça de biopsia - quente
 Biopsitang - varm
 Hotbiopsietang
 Biyopsi Forsepsi – hot
 생검 겸자 - 고온
 Щипцы биопсийные горячие одноразовые

Biopsy Forceps – hot – enteroscope

Pince à biopsie – chaude - entéroscope **Reorder No. 00711295** (350cm)
 Biopsiezange - für Hochfrequenzstrom - Enteroskop
 Pinza per biopsia – a caldo - enteroscopio
 Fórceps para biopsia - caliente – enteroscopia
 Pinça de biopsia - quente - enteroscópio
 Biopsitang - varm - enteroskop
 Hotbiopsietang - enteroscoop
 Biyopsi Forsepsi - hot – enteroskop
 생검 겸자 - 고온 - 내시경
 Щипцы биопсийные горячие для использования с энтероскопом

Biopsy Forceps – pediatric

Pince à biopsie - pédiatrique **Reorder No. 00711261** (160cm)
 Biopsiezange - Pädiatrie **Reorder No. 00711262** (160cm)
 Pinza per biopsia - pediatrica
 Fórceps para biopsia - uso pediátrico
 Pinça de biopsia - pediátrica
 Biopsitang - pædiatrisk
 Biopsietang – kind
 Biyopsi Forsepsi – pediyatrik
 생검 겸자 - 소아
 Щипцы биопсийные педиатрические

INSTRUCTIONS FOR USE

Instructions d'utilisation
 Gebrauchsanweisung
 Istruzioni per l'uso
 Instrucciones de uso
 Instruções de Utilização
 Brugsanvisning
 Gebruiksaanwijzing
 Kullanım talimatları
 사용 설명서
 Инструкция по применению

This product has been manufactured not to include latex.

INTENDED USE:

Disposable biopsy forceps are used to obtain mucosal tissue samples during endoscopic procedures.

WARNINGS AND PRECAUTIONS:

- Endoscopic biopsy should only be performed by persons having adequate training and familiarity with endoscopic biopsy techniques. Consult the medical literature relative to techniques, contraindications, complications and hazards prior to any endoscopic procedure. If reference information for endoscopic technique or contraindications is needed call US Endoscopy.
- Do not use this device if there is any apparent shipping or handling damage.
- Endoscopic biopsy should not be attempted until proficiency has been developed with biopsy forceps.
- Excessive force should not be necessary and should be avoided as not to damage the forceps.
- Biopsy forceps should never be forced through the endoscope's channel. If resistance is met due to excess angulation of the endoscope, it may be necessary to decrease the angulation to allow passage of the forceps.
- Biopsy forceps should not be extracted through the biopsy channel unless the biopsy cups have been completely closed (with or without a specimen.)
- US Endoscopy did not design this device to be reprocessed or reused, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of the device to ensure patient and/or user safety.
- **These disposable medical devices are not intended for reuse.** Any institution, practitioner or third party who reprocesses, refurbishes, remanufactures, resterilizes and/or reuses these disposable medical devices must bear full responsibility for product safety and effectiveness.
- **For product numbers 00711211, 00711295:** Endoscopic hot biopsy should only be performed by persons having a thorough understanding of the principals of diathermic energy. Consult the electro-surgical generator manufacturer's instruction booklet for proper settings and use of the generator.
- In order to ensure that the insulating properties of the device are not compromised, do not exceed the maximum rated peak voltage of 2500 V. Care should be exercised when using hot biopsy forceps to avoid causing unwanted burning or perforation of the wall of the GI tract.

DIRECTIONS FOR USE:

1. Select the forceps compatible with the biopsy channel of the endoscope.
2. Open the sterile pouch. Remove and inspect forceps for obvious damage. Uncoil the forceps and actuate the handle several times to be sure the biopsy cups open and close properly. Should there be evidence of damage, or if the unit is not functioning properly, **DO NOT USE THIS PRODUCT.**
3. After endoscopic visualization of the area for biopsy, place thumb in loop of forceps handle, rest middle and index finger (placing gentle pressure) on the handle's spool edge.
4. With gentle pressure on the handle to keep the cups closed, insert the forceps through the biopsy channel of the endoscope using short strokes until the cups are visualized exiting the endoscope.
5. Manuever the forceps towards the targeted site. Open the cups by activating the handle. Advance the opened forceps over the targeted site and close the forceps. Use only enough pressure to bite the tissue. Over-exertion could cause forceps to become misaligned or to fail.
6. Continue to apply gentle pressure on the handle as the forceps are withdrawn from the channel.
7. Remove and prepare the tissue specimen according to standard technique for histologic evaluation.

DIRECTIONS FOR USE: (FOR PRODUCT NUMBERS 00711211, 00711295 Only)

1. Follow steps 1 – 4 above.
2. Attach the correct active cord to the handle and plug into generator. Set the generator to the pre-determined settings.
3. Manuever the forceps towards the targeted site. Open the cups by activating the handle. Advance the opened forceps over the targeted site and close the forceps. Use only enough pressure to bite the tissue. Over-exertion could cause forceps to become misaligned or to fail.
4. The time and power settings should be determined clinically to control hemostasis.
5. Once the biopsy has been completed and hemostasis has been attained, the procedure can continue to the next step.
6. Continue to apply gentle pressure on the handle as the forceps are withdrawn from the channel.
7. Remove and prepare the tissue specimen according to standard technique for histologic evaluation.

PRODUCT DISPOSAL:

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

ISSUED DATE: June 2013

Warning:

An issued or revision date for these instructions is included for the user's information. In the event two years has elapsed between this date and product use, the user should contact US Endoscopy to determine if additional information is available.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

US Endoscopy, a wholly owned subsidiary of STERIS Corporation.

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Listening...and delivering solutions®, le logo US Endoscopy listening...and delivering solutions, ainsi que toutes les marques portant le signe ® ou ™ sont déposées auprès de l'U.S. Patent and Trademark Office ou sont des marques de commerce détenues par US Endoscopy.
Listening...and delivering solutions®, das US Endoscopy „listening...and delivering solutions“-Design und alle mit ® oder ™ bezeichneten Marken sind eingetragene Marken, die beim U.S. Patent and Trademark Office (US-Patent- und Markenamt) eingetragen wurden, oder Marken, die sich im Besitz von US Endoscopy befinden.
Listening...and delivering solutions®, o design US Endoscopy listening...and delivering solutions e todas as marcas assinaladas com ® ou ™ estão registadas no Departamento de Registos de Patentes dos EUA (U.S. Patent and Trademark Office) ou são marcas comerciais detidas pela US Endoscopy
Listening...and delivering solutions®, US Endoscopy listening...and delivering solutions design og alle varemærker markeret med ® eller ™ er registreret i USA's patentregister (U.S. Patent and Trademark Office) eller er varemærker, der tilhører US Endoscopy.
Listening...and delivering solutions[®], het ontwerp listening...and delivering solutions van US Endoscopy en alle aanduidingen met [®] of [™] zijn geregistreerd bij het Patent and Trademark Office in de Verenigde Staten of zijn handelsmerken van US Endoscopy.
Listening...and delivering solutions®, US Endoscopy listening...and delivering solutions의 디자인과 ® 또는 ™ 이 표시된 모든 마크는 미국 특허청(the U.S. Patent and Trademark Office)에 등록되었거나 US Endoscopy 소유의 등록 상표입니다.