## Roth Net® Platinum™ – food bolus

Roth Net® Platinum™ – bol alimentaire  
Roth Net® Platinum™ – Nahrungsmittelbolus  
Roth Net® Platinum™ – retino boli alimentari  
Roth Net® Platinum™ – extracción de bolos de comida  
Roth Net® Platinum™ – bólus alimentar  
Roth Net® Platinum™ – fødebolus  
Roth Net® Platinum™ – voedselbolus  
Roth Net® Platinum™ – gıda kitlesi  
Roth Net® Platinum™ – 음식 알약

## Roth Net® Platinum™ – universal

Roth Net® Platinum™ – universel  
Roth Net® Platinum™ – Universal  
Roth Net® Platinum™ – universale  
Roth Net® Platinum™ – universal  
Roth Net® Platinum™ – universal  
Roth Net® Platinum™ – universal  
Roth Net® Platinum™ – universeel  
Roth Net® Platinum™ – evrensel  
Roth Net® Platinum™ – 범용

## Roth Net® Platinum™ – polyp

Roth Net® Platinum™ – polype  
Roth Net® Platinum™ – Polyp  
Roth Net® Platinum™ – retino per polipi  
Roth Net® Platinum™ – extirpación de pólipos  
Roth Net® Platinum™ – pólipo  
Roth Net® Platinum™ – polyp  
Roth Net® Platinum™ – poliep  
Roth Net® Platinum™ – polip  
Roth Net® Platinum™ – 용종

### Reorder No. 00711155

Référence de commande 00711155  
Nachbestell-Nr. 00711155  
N. di riordino 00711155  
N.º de pedido 00711155  
N.º de encomenda 00711155  
Genbestillingsnr. 00711155  
Nabestelnr. 00711155  
Yeni Sipariş Numarası: 00711155  
주문 번호00711155

### Reorder No.00715050

Référence de commande 00715050  
Nachbestell-Nr. 00715050  
N. di riordino 00715050  
N.º de pedido 00715050  
N.º de encomenda 00715050  
Genbestillingsnr. 00715050  
Nabestelnr. 00715050  
Yeni Sipariş Numarası: 00715050  
주문 번호00715050

### Reorder No. 00715150

Référence de commande 00715150  
Nachbestell-Nr. 00715150  
N. di riordino 00715150  
N.º de pedido 00715150  
N.º de encomenda 00715150  
Genbestillingsnr. 00715150  
Nabestelnr. 00715150  
Yeni Sipariş Numarası: 00715150  
주문 번호00715150

---

## INSTRUCTIONS FOR USE

Instructions d'utilisation  
Gebrauchsanweisung  
Istruzioni per l'uso  
Instrucciones de uso  
Instruções de Utilização  
Brugsanvisning  
Gebruiksaanwijzing  
Kullanım talimatları  
사용 설명서

---

*listening...and delivering solutions\**

5976 Heisley Road  
Mentor, OH 44060 USA

phone +1 440 / 639.4494

fax +1 440 / 639.4495

email [global@usendoscopy.com](mailto:global@usendoscopy.com)

[www.usendoscopy.com](http://www.usendoscopy.com)

NON-STERILE

Rx Only (U.S.A.)

### Intended Use:

The disposable Roth Net<sup>®</sup> retrieval device product line is used in the endoscopic retrieval of foreign body, food bolus and excised tissue such as polyps.

<u>Description</u>	<u>Product Number</u>	<u>Sheath Diameter</u>	<u>Sheath Length</u>
Roth Net <sup>®</sup> Platinum™ – food bolus	00711155	2.5mm	160cm
Roth Net <sup>®</sup> Platinum™ – universal	00715050	2.5mm	230cm
Roth Net <sup>®</sup> Platinum™ – polyp	00715150	3.0mm	230cm

### Warnings and Precautions:

- Do not use this device if there is any apparent shipping or handling damage.
- The endoscopic retrieval of foreign objects, food bolus, or polyps should only be performed by persons having adequate training and familiarity with endoscopic techniques. Consult the medical literature relative to techniques, complications and hazards prior to the performance of any endoscopic procedure. If reference information for endoscopic techniques utilizing the Roth Net<sup>®</sup> retrieval device product line is desired/needed, contact US Endoscopy; case reports describing techniques are available.
- Avoid blindly passing the Roth Net<sup>®</sup> device past any foreign object, particularly if the entire lumen is blocked.
- Care should be exercised when grasping the targeted object, material, or polyp to be retrieved to avoid inadvertently grasping tissue or organs not intended for retrieval.
- Keep gentle traction on the device during retrieval so that the retrieved object does not become loosened, dislodged or aspirated into the trachea.
- Optimal technique using a Roth Net<sup>®</sup> device to remove a food bolus or any foreign object will vary depending upon the type of bolus/object or clinical conditions faced during retrieval (eg. bolus position, consistency, anatomical, or disease issues such as strictures, bolus-induced ulcerations, etc.). The use of the Roth Net<sup>®</sup> device or any technique described in any document or publication is at the discretion of the medical personnel attempting the retrieval.
- **The Roth Net<sup>®</sup> Platinum™ series is designed using a flat wire. Closing or retracting the handle too tightly may cause damage to the device's flat wire. Do not exert excessive "bending" pressure on the open (deployed) device; doing so may disfigure either the flat wire and/or the wire frame's shape.**
- After extraction, if retrieved material does not exit the net upon net deployment, the material should be removed from the net by rinsing the net in a basin of water; **do not "pick" retrieved material from the net** with a finger as this may damage the net.
- These devices are not recommended for the retrieval of sharp foreign objects.
- The following conditions may not allow the Roth Net<sup>®</sup> device to function properly: (1) advancing the handle to the open position with too much speed or force, (2) attempting to pass or open the device in an extremely articulated endoscope, (3) attempting to actuate the device in an extremely coiled position and/or (4) actuating the device when the handle is at an acute angle in relation to the sheath.
- Short strokes, 1"-1.5" (2.5cm – 3.8cm) in length, are recommended throughout device passage to avoid sheath kinking.
- US Endoscopy did not design this device to be reprocessed or reused, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of the device to ensure patient and/or user safety.
- The Roth Net<sup>®</sup> Platinum™ retrieval device product line is intended for single patient use only. Any institution, practitioner, or third party who reprocesses, refurbishes, remanufactures, reesterilizes, and/or reuses these disposable devices must bear full responsibility for their safety and effectiveness.

### Directions for Use:

1. Prior to clinical use you should familiarize yourself with the device.
  - Read the "Warnings and Precautions" and the "Directions for Use." Review the diagrams.
2. To repack the Roth Net<sup>®</sup> device into the sheath for introduction into the endoscope:
  - Uncoil the entire device and drape in a "U" shaped configuration, holding the proximal device end in one hand, and the distal sheath end in the opposite hand.
  - With the device in this position, carefully retract the handle until the Roth Net<sup>®</sup> device is completely withdrawn into its sheath (see Fig. 3).
3. Check that the accessory channel diameter of the endoscope is compatible with the diameter of the Roth Net<sup>®</sup> sheath.
4. Once the object has been endoscopically identified, advance the retracted Roth Net<sup>®</sup> device into the accessory channel of the endoscope using short strokes until the distal end of the sheath is endoscopically visualized.
5. Advance the distal end of the sheath slightly past the object, bolus, or polyp to be retrieved, however, avoid blindly passing the device past any object, bolus, or polyp if the entire lumen is blocked. Open the Roth Net<sup>®</sup> device by advancing the handle forward until it stops. Confirm the net is fully open via endoscopic observation (see Fig.1). **The following conditions may cause the device to function improperly: (1) advancing the handle to the open position with too much speed or force, (2) attempting to pass or open the device in an extremely articulated endoscope, (3) attempting to actuate the device in an extremely coiled position and/or (4) actuating the device when the handle is at an acute angle in relation to the sheath.**
6. Endoscopically manipulate the Roth Net<sup>®</sup> device over the object, bolus, or polyp. In most instances, the object should be centered and proximal in the net (for best results) prior to net closure, **however for food bolus retrieval techniques may vary.**
7. **Food bolus retrieval**
  - Optimal technique using a Roth Net<sup>®</sup> device to remove a food bolus or any foreign object will vary depending upon the type of object/bolus or clinical conditions faced during retrieval (eg. bolus position, consistency, anatomical, or disease issues such as strictures, bolus-induced ulcerations, etc.).
  - One observed technique demonstrated the clinician sufficiently insufflating the esophagus, then positioning the tip of the device sheath alongside the *outer edge* of the food bolus while the assistant partially opened the net; the clinician then simultaneously, but gently *withdrew the proximal end of the device .5 to 1.5 cm from the biopsy inlet* (this maneuver resulted in positioning the open net midline to the bolus). The tip of the endoscope was then deflected to drape the net across the proximal face of the bolus. Once this was accomplished, the net was gently closed until *mild* resistance was felt. Removal of the bolus was endoscopically observed, confirming sufficient capture of the bolus throughout its extraction. If and when repositioning of the net on the bolus was necessary (distal to the tracheal bifurcation), the net was fully opened and the capture sequence described above was repeated.
  - **Closing or retracting the handle too tightly may cause damage to the device's flat wire. Do not exert excessive "bending" pressure on the open (deployed) device; doing so may disfigure either the flat wire and/or the wire frame's shape.**
  - Retrieved food should be removed from the net by rinsing the net in a basin of water; do not "pick" retrieved food from the net with a finger as this may damage the net.

8. Once entrapment is accomplished, retract the handle until it stops (see Fig. 2). Continuous gentle traction should be applied on the handle to keep the Roth Net<sup>®</sup> device closed. Excess pressure may cause the net to rupture or tear. Retract the device, retrograde, toward the endoscope so the Roth Net<sup>®</sup> device is endoscopically observed yet not obscuring the endoscopic view. Endoscopic observation is necessary during extubation as not to lose sight of the object, bolus, or polyp during removal.
9. Once the endoscope and the Roth Net<sup>®</sup> device have been extubated, the retrieved object, bolus, or polyp can be removed from the net by advancing the handle forward to open the Roth Net<sup>®</sup> device (see Fig. 1). Specimens should be prepared according to standard technique for histologic evaluation.
10. **To repack the Roth Net<sup>®</sup> device into the sheath for introduction to or removal from the endoscope or for additional retrievals:**
  - Retract the handle until the Roth Net<sup>®</sup> device is completely withdrawn into its sheath (see Fig. 3).
11. Once the endoscopic retrieval has been completed, remove the device from the endoscope.

**Product Disposal:**

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.




**Issued Date: April 2010**

**Warning:**

An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact US Endoscopy to determine if additional information is available.

**Listing...and delivering solutions<sup>®</sup>**, the US Endoscopy listing...and delivering solutions design, and all marks denoted with <sup>®</sup> or <sup>™</sup> are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy, U.S. patent numbers 6,814,739 and 5,643,283, other patents pending

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Listening...and delivering solutions®, o design US Endoscopy listening...and delivering solutions e todas as marcas assinaladas com ® ou ™ estão registadas no Departamento de Registos de Patentes dos EUA (U.S. Patent and Trademark Office) ou são marcas comerciais detidas pela US Endoscopy, números de patentes nos E.U.A. 6,814,739 e 5,643,283, outras patentes pendentes
**Listening...and delivering solutions®**, US Endoscopy listening...and delivering solutions 디자인, 그리고® 또는 ™ 으로 표시된 모든 상표는 미국 특허청 (U.S. Patent and Trademark Office) 에 등록되어 있거나 US Endoscopy가 소유하고 있는 상표입니다. 미국 특허 번호 6,814,739 및 5,643,283, 기타 특허 출원 중