### Articulator™ - injection needle

Articulator™ - aiguille pour injection  
Articulator™ - Injektionsnadel  
Articulator™ - ago per iniezione  
Articulator™ - aguja de inyección  
Articulator™ - agulha de injeção  
Articulator™ - injektionsnål  
Articulator™ - injectienaald  
Articulator™ - enjeksiyon iğnesi  
Articulator™ - 주사기  
Articulator™ – инъекционная игла

Reorder No. 00711803  
Reorder No. 00711804  
Reorder No. 00711807  
Reorder No. 00711808  
Reorder No. 00711810

(4mm / 160cm / 25ga / 2.5mmØ)  
(5mm / 160cm / 25ga / 2.5mmØ)  
(5mm / 230cm / 25ga / 2.5mmØ)  
(5mm / 350cm / 25ga / 2.5mmØ)  
(5mm / 230cm / 25ga / 2.5mmØ)

### Carr-Locke - injection needle

Carr-Locke - aiguille pour injection  
Carr-Locke - Injektionsnadel  
Carr-Locke - ago per iniezione  
Carr-Locke - aguja de inyección  
Carr-Locke - agulha de injeção  
Carr-Locke - injektionsnål  
Carr-Locke - injectienaald  
Carr-Locke - enjeksiyon iğnesi  
Carr-Locke – 주사기  
Carr-Locke – инъекционная игла

Reorder No. 00711811  
Reorder No. 00711812  
Reorder No. 00711813  
Reorder No. 00711814  
Reorder No. 00711822  
Reorder No. 00711823  
Reorder No. 00711824

(5mm / 230cm / 25ga / 2.5mmØ)  
(5mm / 230cm / 23ga / 2.5mmØ)  
(5mm / 230cm / 25ga / 2.5mmØ)  
(4mm / 230cm / 23ga / 2.5mmØ)  
(5mm / 230cm / 25ga / 1.8mmØ)  
(5mm / 230cm / 25ga / 2.5mmØ)  
(5mm / 230cm / 23ga / 2.5mmØ)

### Vari-Safe™ - injection needle

Vari-Safe™ - aiguille pour injection  
Vari-Safe™ - Injektionsnadel  
Vari-Safe™ - ago per iniezione  
Vari-Safe™ - aguja de inyección  
Vari-Safe™ - agulha de injeção  
Vari-Safe™ - injektionsnål  
Vari-Safe™ - injectienaald  
Vari-Safe™ - enjeksiyon iğnesi  
Vari-Safe™ - 주사기  
Vari-Safe™ – инъекционная игла

Reorder No. 00711819

(5mm / 230cm / 23ga / 2.5mmØ)

---

## INSTRUCTIONS FOR USE

MODE D'EMPLOI  
GEBRAUCHSANLEITUNG  
ISTRUZIONI PER L'USO  
INSTRUCCIONES DE USO  
INSTRUÇÕES DE UTILIZAÇÃO  
BRUGSANVISNING  
GEBRUIKSAANWIJZING  
KULLANIM TALİMATLARI  
사용 설명서  
ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ

---

*listening... and delivering solutions\**

5976 Heisley Road  
Mentor, OH 44060 USA

phone +1 440 / 639.4494

fax +1 440 / 639.4495

email [global@usendoscopy.com](mailto:global@usendoscopy.com)

[www.usendoscopy.com](http://www.usendoscopy.com)

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***This product has been manufactured not to include latex.***

#### **Intended Use:**

These single use, disposable, flexible sheath injection needles are intended to be used for the injection of various types of media through flexible endoscopes.

#### **Warnings and Precautions:**

- Flexible endoscopy procedures should only be performed by persons having adequate training and familiarity with endoscopic techniques. Consult the medical literature relative to techniques, contraindications, complications and hazards prior to any endoscopic procedure.
- These injection needles are single use, disposable medical devices; they are not intended for reuse. Any institution, practitioner or third party who reprocesses, refurbishes, remanufactures, sterilizes, and/or reuses these disposable devices must bear full responsibility for their safety and effectiveness.
- US Endoscopy did not design these devices to be reprocessed or reused, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of the device to ensure patient and/or user safety.
- **For Carr-Locke injection needles only:** Do not attempt to actuate the injection needle with the catheter in a straightened position. The injection needle has been engineered and manufactured with a specific "preload" that is necessary to assure the full projection and full retraction of the needle from the distal hub when the device is in very tortuous configurations. If the injection needle is actuated with the catheter in a straightened position, the "preload" may cause the needle's inner catheter to collapse and not allow for adequate injection. When handling the injection needle, always be aware of the distal tip of the needle, and visually confirm that the needle is fully retracted into the device.
- **For Vari-Safe™ injection needles only:** Do not attempt to activate or test the injection needle outside of the endoscope if the device is coiled more than once. It is possible that this may cause the injection needle to become lodged proximal to the distal tip spring.
- **Compatibility with media:**
  - All of the injection needles included in this IFU (Instructions for Use) are compatible with the following medicinal substances and solutions: Sodium Morrhuate, Sodium Tetradecyl Sulfate, Ethynolamine Oleate, Epinephrine, Botulinum Toxin Type A, Polidocanol, Ethanol, Hypertonic Saline and Staining Media. **Note: Sodium Morrhuate, Sodium Tetradecyl Sulfate, Polidocanol, and Ethanol have not been approved by the United States Federal Food and Drug Administration and should not be used in endoscopic procedures.**
  - The Carr-Locke 23 gauge injection needle (product # 00711812, #00711814, #00711824) and the Vari-Safe™ injection needle (product # 00711819) are compatible with Triamcinolone.
  - The Carr-Locke 23 gauge injection needle (product # 00711812, #00711814, #00711824) is compatible with Hyaluronic Acid, Hydroxypropyl Methylcellulose, Glycerine, Lipiodol, and N-butylcyanocrylate. **Note: Lipiodol and N-butylcyanocrylate have not been approved by the United States Federal Food and Drug Administration and should not be used in endoscopic procedures.**
  - The compatibility of the injection needles with other media has not been verified and therefore should not be attempted.
- **Device Compatibility: The Carr-Locke injection needle (product number 00711813) is the only needle included in the IFU that is compatible with the Pentax Model 90 series High Definition gastroscopes.**

**Contraindications:** Contraindications include those specific to any endoscopic procedure, as well as those specific to endoscopic needle therapy.

#### **Directions for Use:**

1. (a) **For Articulator™ & Carr-Locke Injection Needles Only:** Carefully remove the tip protector from the distal tip of the needle. Actuate the proximal luer back and forth, visually observing the needle being adequately advanced from the device's distal end. The spring-loaded handle assists in needle retraction; however, always visually check the retraction as well. **(Note: Product #00711808 version used for enteroscopy does not have this spring-loaded handle feature).**  
(b) **For Carr-Locke Injection Needle Usage in a Duodenoscope:** When using this needle in a duodenoscope, before actuating the elevator to the full position, assure that the needle sheath is projected adequately past the end of the scope to prevent damage to the needle which, in the retracted position, is positioned inside the spring, proximal to the distal tip.  
(c) **For Vari Safe™ injection needles only:** Carefully remove the tip protector from the distal tip of the needle. Uncoil the needle sheath and ensure that the needle is properly projected. Unlock the luer lock fitting and actuate the proximal luer back and forth, visually observing the needle being adequately advanced from the distal end. The return spring assists in the needle retraction; however, always visually check the retraction as well.  
**Note: Do not stretch the device because stretching may cause the needle projection length to be altered.**
2. To lock the needle in the out position, gently press the proximal luer until it meets the handle's luer, and then slowly twist the proximal luer clockwise until the two luers are fully engaged.
3. Using accepted medical methods, purge the needle assembly of air using 0.5-0.75cc of injection media.
4. Ensure that the needle is fully retracted prior to inserting it into the endoscope. To unlock the needle, slowly twist the proximal luer counterclockwise until the two luers disengage.
5. Proceed with injection of appropriate media using approved medical techniques until desired results are accomplished.
6. When moving from one injection site to the next, it is advisable to retract the needle into the sheath to protect the mucosa, as well as the endoscope, from trauma/damage.
7. When removing the injection needle from the endoscope, fully retract the needle into the sheath to protect the patient and clinician from injury, and to protect the endoscope from potential damage.

#### **Product Disposal:**

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

**Issued Date:** January 2017

#### **Warning:**

An issued or revision date for these instructions is included for the user's information. In the event two years has elapsed between this date and product use, the user should contact US Endoscopy to determine if additional information is available.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

This product and/or its use is protected under United States patent number 6,666,847.

US Endoscopy, a wholly owned subsidiary of STERIS Corporation.

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Listening...and delivering solutions®, le logo US Endoscopy listening...and delivering solutions, ainsi que toutes les marques portant le signe ® ou ™ sont déposées auprès de l'U.S. Patent and Trademark Office ou sont des marques commerciales d'US Endoscopy.
Listening...and delivering solutions®, das US Endoscopy listening...and delivering solutions-Design und alle mit ® oder TM bezeichneten Marken sind beim U.S. Patent und Trademark Office (US-Patent- und Markenamt) eingetragene Marken oder Warenzeichen, die sich im Besitz von US Endoscopy befinden.
Listening...and delivering solutions®, o design US Endoscopy listening...and delivering solutions e todas as marcas assinaladas com ® ou ™ estão registadas no Departamento de Registos de Patentes dos EUA (U.S. Patent and Trademark Office) ou são marcas comerciais detidas pela US Endoscopy.
Listening...and delivering solutions®, US Endoscopy listening...and delivering solutions design og alle varemærker markeret med ® eller ™ er registreret i USA's patentregister (U.S. Patent and Trademark Office) eller er varemærker, der tilhører US Endoscopy.
Listening...and delivering solutions®, het ontwerp listening...and delivering solutions van US Endoscopy en alle aanduidingen met ® of ™ zijn geregistreerd bij het Patent and Trademark Office in de Verenigde Staten of zijn handelsmerken van US Endoscopy.
Listening...and delivering solutions®, US Endoscopy listening...and delivering solutions 디자인, 그리고 ® 또는 ™으로 표시된 모든 상표는 미국 특허청(U.S. Patent and Trademark Office)에 등록되어 있거나 US Endoscopy가 소유하고 있는 상표입니다.