Histoguide™ wire-guided forceps

Pince guide par fil Histoguide™

Histoguide™ drahtgeführte Zange

Pinze filo-guidate Histoguide™

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*listening...and delivering solutions**

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639.4494

fax +1 440 / 639.4495

email global@usendoscopy.com

www.usendoscopy.com

Not made with natural rubber latex.

Intended Use: The disposable wire-guided forceps is intended to sample tissue from lesions that can occur within the gastrointestinal tract.

Description	Product Number	Sheath Diameter	Jaw Diameter	Jaw Opening Width	Device Length	Combined Wire Guide & Sheath Diameter (See Fig. 2)	Recommended Guidewire Diameter
Histoguide™ wire-guided forceps	00711660	2.1mm	2.2 mm	8 mm	2300 mm	0.140in / 3.5mm	0.035in / 0.89mm

Design Features:

- **Easy to use spring loaded handle is designed to keep the jaws closed until ready for use and prevent specimen loss during exchanges. SEE FIG.1.**
- **Wire guide for easy and safe guide wire exchange. SEE FIG.2.**
- **Metallic biopsy cups for enhanced fluoroscopic imaging.**

Warnings and Precautions:

- Consult the medical literature relative to techniques, complications and hazards prior to the performance of any endoscopic procedure and ERCP.
- The sampling of tissue within and outside the gastrointestinal tract should only be performed by persons having adequate training and familiarity with endoscopic techniques procedures and ERCP.
- Care should be exercised when grasping the targeted tissue to avoid inadvertently grasping tissue or organs not intended for retrieval.
- For best results keep the guide wire wet.
- If any resistance is felt during device manipulation **STOP** and determine the cause prior to proceeding.
- Fluoroscopic guidance is recommended when using this device.
- Sphincterotomy is recommended when using this device.
- Optimal technique for tissue sampling will vary depending upon the type or clinical conditions found during endoscopy and ERCP.
- Use caution to minimize risk of any SHARPS injury.
- Do not use in the pancreas due to pancreatitis risk.
- Wear appropriate Personal Protective Equipment.
- Do not attempt to reuse, reprocess, refurbish, remanufacture or resterilize this device. US Endoscopy did not design this device nor is it intended to be reused, reprocessed, refurbished, remanufactured, or resterilized. Performing such activities on this disposable medical device presents a safety risk to patients (i.e. compromised device integrity, cross-contamination, infection).

Complications:

- Those associated with gastrointestinal endoscopy include but are not limited to:
 - Pancreatitis
 - Perforation.
 - Bleeding or hemorrhage.
 - Aspiration, fever, infection, allergic reaction to medication.
 - Hypotension, respiratory depression or arrest cardiac arrhythmia or arrest.

Contraindications:

- Those specific to any endoscopic procedure.
- Those specific to ERCP procedure or procedure where gaining access to a specific site is desired.
- Patients with known or suspected coagulopathy.
- Use in the pancreas

Prior to use:

1. Confirm that the biopsy forceps are compatible with the inner channel of the endoscope (**min 4.2mm channel**).
2. Read the entire Instructions for Use and familiarize yourself with the device.
3. Visually inspect the packaging for damage.
4. Open the package, remove and uncoil the entire device slowly.
5. Visually inspect the device for damage.
6. Drape the sheath in a “U” shaped configuration.
7. Actuate the device by moving the Slider back and forth to confirm that the jaws open and close smoothly. **SEE FIG.1. If the unit does not function properly, or there is evidence of damage to the device or package (e.g. bends, kinks, or misaligned jaws, damaged packaging) do not use this product and contact your local Product Specialist.**

Directions for Use:

NOTE: Once access to the sampling site is achieved and a guidewire is placed through the stricture, proceed with the following steps;

1. Identify the lesion to be sampled via fluoroscopic imaging, and **maintain visualization throughout the use of the forceps.**
2. Insert the guidewire into the WIRE GUIDE attached to the biopsy cups. **SEE FIG.2.**
3. Slide the forceps over the guidewire.

- Squeeze the Slider and Thumb Ring to prevent inadvertent opening of the biopsy cups when inserting the instrument into the endoscope.

Note: Spring loaded handle is designed to minimize biopsy cups from opening during the exchange process. **SEE FIG.1.**

- Advance the forceps into the endoscope biopsy channel using short strokes (1.0 to 1.5 inches) to avoid kinking the sheath.

Caution: Do not force the forceps through the endoscope. Forcing the device through the endoscope could cause damage to the device or endoscope.

- Maintain fluoroscopic imaging.
- Advance the forceps out of the endoscope and advance the device to desired area within the lesion.

Note: Reduce the angulation of the endoscope or lower the forceps elevator in the endoscope if resistance is felt.

- OPEN the cups by moving the Slider away from the Thumb Ring. **SEE FIG.1.**
- ADVANCE the forceps with the cups OPEN to the desired tissue while maintaining fluoroscopic imaging.
- CLOSE the cups by moving the Slider toward the Thumb Ring, and grasp the tissue. **SEE FIG.1.**

Caution: Use only enough pressure to bite the tissue. Over-exertion could cause forceps to become misaligned or to fail.

- Maintain pressure on the Slider and Thumb Ring.
- Grasp the biopsy forceps metal sheath.
- Pull in a direction to remove the tissue sample from the lesion.
- Remove the forceps from endoscope using slow continuous strokes to minimize the scatter of bodily fluids.
- Remove the forceps from the guidewire.
- Transfer specimen to a specimen jar and medium for preservation for pathology per institutional guidelines.
- Rinse the forceps in sterile water while actuating handle until any excess tissue/specimen is removed.
- If additional tissue sample is desired, visually inspect the device for damage and then repeat steps 1-15.


FIG. 1


FIG. 2

Issued Date: May 2018

Product Disposal:

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

WARNING:

An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact US Endoscopy to determine if additional information is available.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

US Endoscopy, a wholly owned subsidiary of STERIS Corporation

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Listening...and delivering solutions®, le logo US Endoscopy listening...and delivering solutions, ainsi que toutes les marques portant le signe ® ou ™ sont des marques déposées auprès de l'U.S. Patent and Trademark Office ou sont des marques commerciales détenues par US Endoscopy.
Listening...and delivering solutions®, das „listening...and delivering solutions“-Design und alle mit ® oder ™ versehenen Marken sind beim United States Patent and Trademark Office (Patentamt der Vereinigten Staaten) registriert oder befinden sich im Besitz von US Endoscopy
Listening...and delivering solutions®, o design US Endoscopy listening...and delivering solutions e todas as marcas assinaladas com ® ou ™ estão registadas no Departamento de Registos de Patentes dos EUA (U.S. Patent and Trademark Office) ou são marcas comerciais detidas pela US Endoscopy.
Listening...and delivering solutions®, designet af US Endoscopy listening...and delivering solutions, og alle mærker markeret med ® eller ™ er registreret hos U.S. Patent and Trademark Office eller varemærker ejet af US Endoscopy.
Listening...and delivering solutions®, het ontwerp Listening...and delivering solutions van US Endoscopy en alle aanduidingen met ® of ™ zijn gedeponeerd bij het Patent and Trademark Office in de Verenigde Staten of zijn handelsmerken van US Endoscopy
Listening...and delivering solutions®, US Endoscopy listening...and delivering solutions 디자인, 그리고 ® 또는 ™으로 표시된 모든 상표는 미국 특허청(U.S. Patent and Trademark Office)에 등록되어 있거나 USEndoscopy가 소유하고 있는 상표입니다.