iSnare[®] system – Lariat[®] snare
Système iSnare[®] – anse Lariat[®]
iSnare[®] System – Lariat[®] Schlinge
Sistema iSnare[®] – Ansa Lariat[®]
Asa iSnare system[®] - Lariat[®]
Sistema iSnare[®] – ansa Lariat[®]
iSnare[®]-systemet – Lariat[®]-slyngje
iSnare-systeme[®] - Lariat[®]-poliepsnoerder
iSnare[®] sistemi – Lariat[®] snare
iSnare[®] system - Lariat[®] 스네어
Система iSnare[®] – петля Lariat[®]

Reorder No. 00711084
Référence de commande 00711084
Nachbestell-Nr. 00711084
N. di riordino 00711084
Nº de pedido 00711084
N.º de encomend 00711084
Genbestillingsnr. 00711084
Nabestelnr. 00711084
Yeni Sipariş Numarası 00711084
주문 번호 00711084
№ для заказа 00711084

INSTRUCTIONS FOR USE

MODE D'EMPLOI

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BRUGSANVISNING

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KULLANIM TALIMATLARI

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www.usendoscopy.com

CE 0086

This device is not made with natural rubber latex.

Intended use: The US Endoscopy iSnare® system – Lariat® snare is intended to be used to inject various types of media and to grasp, dissect and transect tissue during gastrointestinal endoscopy procedures. The device can be used with or without monopolar diathermic energy.

Device description:

The iSnare® system – Lariat® snare is a single catheter device incorporating both an injection needle and polypectomy snare to be used through flexible endoscopes. The needle is intended to be used for the injection of various types of media. The snare is intended to grasp, dissect, and transect tissue during gastrointestinal endoscopy procedures. The snare can be used with or without monopolar diathermic energy (hot snaring and cold snaring).

Description	Product Number	Sheath Diameter	Sheath Length	Approximate Snare Size (Width x Length)	Snare Wire Diameter	Needle Size	Active Cord Style
iSnare system® - Lariat® snare	BX00711084	2.5 mm	230 cm	Oval shape 30 x 55mm Hex shape 10 x 28mm Diamond shape 6 x 15mm	0.41mm 0.016inches	25 ga	Olympus

US Endoscopy recommends the use of an electrical surgical generator with the following parameters:

Frequency: 330-850 KHz
Peak Voltage: 4300Vp
Connection: Olympus style active cord connector
Pulse shape:

- Cut mode: unmodulated sinusoidal alternating voltage
- Coagulation mode: pulse- modulated alternating current

Please consult the instructions for use of the electrical surgical generator or contact the manufacturer for more detailed information.

Warnings and Precautions:

- Consult the medical literature relative to techniques, complications and hazards prior to the performance of any endoscopic procedure.
- Endoscopic procedures should only be performed by persons having adequate training and familiarity with endoscopic techniques in snare polypectomy and injection with an endoscopic injection needle.
- A thorough understanding of the technical principles, clinical applications, and risks associated with monopolar diathermic energy and cold snaring (non-electrical) polypectomy and tissue resection is necessary before using this product.
- Endoscopic polypectomy with monopolar diathermic snares should not be performed without a thorough understanding of the principles of diathermic energy.
- US Endoscopy did not design this device to be reprocessed or reused, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of the device to ensure patient and/or user safety.
- This disposable medical device is not intended for reuse. Any institution, practitioner, or third party, who reprocesses, refurbishes, remanufactures, reesterilizes, and/or reuses this disposable medical device must bear full responsibility for its safety and effectiveness.
- These devices are compatible with an endoscope accessory channel of 2.8mm or larger.
- Short strokes, 1"-1.5" (2.5cm – 3.8cm) in length, are recommended throughout device passage to avoid sheath kinking.
- This device is not designed for use in side viewing endoscopes.
- The following conditions may not allow the iSnare® system - Lariat® snare device to function properly:
 - 1) advancing the handle to the open position with too much speed or force,
 - 2) attempting to pass or open the device in an extremely articulated endoscope,
 - 3) attempting to actuate the device in an extremely coiled position and/or
 - 4) actuating the device when the handle is at an acute angle in relation to the sheath.
- If there is lack of visible cautery effect, consider the following:
 - 1) Check the active cord for secure connection to the snare handle and the electrosurgical generator.
 - 2) Ensure that the return electrode is secure and properly connected to the patient and the electrosurgical generator.
 - 3) Consult the electrosurgical generator manufacturer's instructions for use for proper settings and use of the generator.
 - 4) **In order to ensure that the insulating properties of the monopolar diathermic snares are not compromised, do not exceed the maximum rated peak voltage of 2500 for cut mode and 2500 for coagulation mode.**

Caution - Output power setting of the electrosurgical unit should be set as low as possible to achieve the intended purpose. If the power setting of the generator is not known, set the power level at the minimum value of the expected range and cautiously increase power to achieve the desired tissue effect.

- Observe for whitening of the tissue as visible evidence of electrocoagulation.
- User should apply caution during diathermic procedures to avoid cutting tissue before the application of cautery.
- To avoid patient injury, do not apply excessive force to the snare while resecting tissue.
- Always be aware of the device's distal tip and visually confirm that the injection needle is fully retracted into the catheter to minimize the risk of accidental injury to the patient or clinician.
- The injection needle should be fully retracted into the catheter during use of the snare to prevent inadvertent injury to the mucosa.

- When moving from one injection site to the next, retract the needle into the sheath to protect the mucosa, as well as the endoscope from trauma/damage.
- The snare should be fully retracted into the catheter during use of the needle to prevent injury to the mucosa.
- Disengage the snare from the polyp if there is a risk of complication.
- Use of a smaller syringe (3cc or 5cc) will assist in an easier injection of a thicker media.

Contraindications:

- Contraindications include, but are not limited to, those specific to any endoscopic procedure.
- Coagulopathy
- Poor bowel prep
- Uncooperative patient
- Consult the appropriate medical personnel regarding use of this device in patients with an implanted electronic device.

Prior to Use:

1. Prior to clinical use, inspect and familiarize yourself with the device. If there is evidence of damage (i.e. deformed or bent snare or needle, inoperable handle, damaged packaging), **do not use this product and contact your local Product Specialist.**
2. Ensure that the endoscope has a channel of 2.8mm or larger before insertion.
3. Remove the device from the package and uncoil the entire device and drape in a "U" shaped configuration, holding the proximal end in one hand and the distal sheath in the opposite hand.
4. **To prepare the injection needle for use:**
 - To lock the needle in the out position, gently advance the proximal luer (A) until it meets the handle's luer (B) and then slowly twist the proximal luer clockwise until the two luers are fully engaged (Figure 1).
 - Attach the syringe filled with the desired injection media to the proximal luer connection. (A) (Syringe not provided)
 - The injection needle should be flushed with 0.5-0.75 cc injection media prior to injection into the tissue. This purges the needle of air.
 - To unlock the needle to the retracted position, slowly twist the proximal luer counterclockwise until the two luers are fully disengaged (Figure 2).
 - Ensure that the needle is fully retracted prior to insertion into the endoscope channel.
 - Ensure that the syringe remains attached to the proximal luer during usage.
5. **To prepare the polypectomy snare for use:**
 - Move the finger rings back and forth to confirm that the snare loop opens and closes smoothly (Figure 3 and 4).
 - Ensure the snare loop is fully retracted into the sheath prior to insertion into the endoscope (Figure 4).
 - Position the patient return electrode and connect to electrosurgical unit following instructions from the electrosurgical unit manufacturer.
 - Attach the respective Olympus style active cord to the snare's diathermic handle connection just prior to excision of polyps or tissue.
 - Consult the electrosurgical generator manufacturer's instructions for use for proper settings and use of the electrical surgical generator.

Directions for Use:

Injection Needle Use:

1. Ensure that the needle and polypectomy snare are retracted prior to insertion into the endoscope channel.
2. After endoscopic identification of the desired injection site, pass the device through the accessory channel of the endoscope.
3. Short strokes, 1"-1.5" (2.5cm – 3.8cm) in length, are recommended throughout device passage to avoid sheath kinking.
4. Once the device is at the injection site, gently advance the proximal luer (A) until it meets the handle's luer (B). To lock the needle in place, slowly twist the proximal luer clockwise until the two luers are fully engaged (Figure 1).
5. Using the attached syringe with the desired media, proceed with injection using approved medical techniques until desired results are accomplished.
6. When moving from one injection site to the next, retract the needle into the sheath to protect the mucosa, as well as the endoscope from trauma/damage.
7. When injection is completed and prior to use of the snare, retract the needle fully into the catheter by twisting the proximal luer counterclockwise until the two luers are separated (Figure 2). The needle is spring loaded and once the two luers are separated, the needle will spring back into the sheath.


FIG. 1


FIG. 2

Hot Snare Use:

A thorough understanding of the technical principles, clinical applications, and risks associated **with the use of monopolar diathermic energy (hot snaring)** for the resection of polyp tissue is necessary before using this product.

1. Position patient return electrode and connect to electrosurgical unit following instructions from the electrosurgical unit manufacturer.
2. Attach the respective Olympus style active cord to the snare's diathermic handle connection prior to excision of polyps or tissue.
3. Consult the electrosurgical generator manufacturer's instructions for use for proper settings and use of the electrosurgical generator.

Caution: In order to ensure that the insulating properties of the monopolar diathermic snares are not compromised, do not exceed the maximum rated peak voltage of 2500 for cut mode and 2500 for coagulation mode.

4. When the polyp or tissue has been endoscopically visualized, deploy the snare by moving the finger rings away from the thumb ring (Figure 3).
5. Confirm that the snare loop is fully open via endoscopic observation (Figure 3).

NOTE: If it is observed that the snare tail is protruding from the catheter (see Figure 5a), slightly move the finger rings towards the thumb ring until the snare tail is inside of the catheter (see Figure 5b).

6. Place the snare around the polyp or tissue being resected and resect it using proper endoscopic technique.
7. Once the polyp(s) or tissue has been resected, the snare can be retracted into the sheath by moving the finger rings towards the thumb ring (Figure 4). The resected polyp(s) or tissue should be removed and prepared according to standard technique for histologic evaluation.

Cold Snare Use:

A thorough understanding of the technical principles, clinical applications, and risks associated **without the use of monopolar diathermic energy (cold snaring)** for the resection of polyp tissue is necessary before using this product.

1. If it is intended to resect polyps or tissue *without* the use of monopolar diathermic energy (cold snaring), the use of the active cord and placement of the return electrode is not necessary. If it is then desired to resect *with* the use of monopolar diathermic energy, follow Hot Snare Use.
2. Snaring without the use of monopolar diathermic energy is recommended for the resection of diminutive polyp tissue.
3. When the polyp or tissue has been endoscopically visualized, deploy the snare by moving the finger rings away from the thumb ring (Figure 3).
4. Confirm that the snare loop is fully open via endoscopic observation (Figure 3).

NOTE: If it is observed that the snare tail is protruding from the catheter (see Figure 5a), slightly move the finger rings towards the thumb ring until the snare tail is inside of the catheter (see Figure 5b).

5. Place the snare around the polyp or tissue being resected and resect it using proper endoscopic technique.
6. Once the polyp(s) or tissue has been resected, the snare can be retracted into the sheath by moving the finger rings towards the thumb ring (Figure 4). The resected polyp(s) or tissue should be removed and prepared according to standard technique for histologic evaluation.

FIG. 3


FIG. 4


FIG. 5A

FIG. 5B

Product Disposal:

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Issued Date: October 2016

Warning:

An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact US Endoscopy to determine if additional information is available.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

US Endoscopy, a wholly owned subsidiary of STERIS Corporation.

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Listening...and delivering solutions®, le logo US Endoscopy listening...and delivering solutions, ainsi que toutes les marques portant le signe ® ou ™ sont des marques déposées auprès de l'U.S. Patent and Trademark Office ou sont des marques commerciales détenues par US Endoscopy.
Listening...and delivering solutions®, das US Endoscopy listening...and delivering solutions-Design und alle mit ® oder ™ bezeichneten Marken sind beim U.S. Patent and Trademark Office (US-Patent- und Markenamt) eingetragen, oder Marken, die sich im Besitz von US Endoscopy befinden.
Listening...and delivering solutions®, o design US Endoscopy listening...and delivering solutions e todas as marcas assinaladas com ® ou ™ estão registadas no Departamento de Registos de Patentes dos EUA (U.S. Patent and Trademark Office) ou são marcas comerciais detidas pela US Endoscopy.
Listening...and delivering solutions[®], US Endoscopy listening...and delivering solutions design og alle varemærker markeret med [®] eller [™] er registreret i USA's patentregister (U.S. Patent and Trademark Office) eller er varemærker, der tilhører US Endoscopy.
Listening...and delivering solutions[®], het ontwerp Listening...and delivering solutions van US Endoscopy en alle aanduidingen met [®] of [™] zijn gedeponeerd bij het Patent and Trademark Office in de Verenigde Staten of zijn handelsmerken van US Endoscopy.
Listening...and delivering solutions®, US Endoscopy listening...and delivering solutions 디자인, 그리고 ® 또는 ™으로 표시된 모든 상표는 미국 특허청(U.S. Patent and Trademark Office)에 등록되어 있거나 US Endoscopy가 소유하고 있는 상표입니다.