





US Endoscopy polypectomy portfolio

reference guide

Snare/Device	Exacto® cold snare BX00711115	Lariat® snare BX00711119	Histolock® resection device BX00711117	iSnare® system - Lariat® snare BX00711084
Polyp Type	 Diminutive (1-9mm)	 Multiple Polyps (all types)	 Very Flat/Sessile	 Large Flat/Sessile
Snare Width	9mm	6mm, 10mm, 30mm	14mm	6mm, 10mm, 30mm
Snare Wire Diameter	Braided .012in. (.30mm) 33% thinner than standard snares (.018") for a clean cut	Braided .016in. (.41mm)	Monofilament flat wire	Braided .016in. (.41mm)
Features	<ul style="list-style-type: none"> • 90% of polypectomies are less than 9mm in size¹. • A dedicated cold snare that has a thinner and stiffer wire may engage and hold healthy polyp margin more easily than a standard snare². • Size and shape offer more control for placement around small polyps. 	<ul style="list-style-type: none"> • Three distinct shapes in one device. • Offers versatility to handle cases with large polyps and multiple polyps of varying types and sizes. • Can be used for hot and cold polypectomy. 	<ul style="list-style-type: none"> • Monofilament stiff wire easily recruits and resects difficult flat polyps. • Ideal for challenging EMR cases with scarred or adherent polypoid tissue. • Grabs tissue predictably for a precise resection. 	<ul style="list-style-type: none"> • Injection needle and snare in one device operated by a multi-functional handle. • Ideal for saline assisted polypectomy cases. • Three distinct shapes in one device.

¹Tolliver KA, Rex DK. Colonoscopic Polypectomy. Gastroenterol Clin N Am 37. 2008; 229-251.

²Daniel von Renteln, MD and Heiko Pohl, MD. Pushing the Limit: How to Get the Most Out of Cold Snares. American Journal of Gastroenterology; 111. July 2016; 1217-1219.

Refer to the current **Instructions For Use** for complete product usage instructions.

To order, please call US Endoscopy Customer Service at 1-800-769-8226.

US Endoscopy is a subsidiary of STERIS Corporation.
 © 2017 US Endoscopy or its affiliates. All rights reserved.
 Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.