

Exacto[®]
cold snare




Clinically proven, purpose-built cold snare

Clinically proven to achieve a significantly higher rate of complete resection, the unique design of the Exacto® cold snare offers more control and placement for a precise, clean cut. ^{1,2}

Specially designed to optimize the cold snare technique, the Exacto® cold snare...

- Enables a clean cut with a thin wire design (.012"/0.30mm) while maintaining its shape and stiffness
- Features a unique, shield shape that maximizes its width for optimal control and placement
- Aids in safe and effective resection with a 9mm mini snare size that reduces the risk of creating large defects
- Offers a thin, stiff catheter that decreases buckling while aiding in placement



For more information on the Exacto® cold snare, contact US Endoscopy today.

product number	description	sheath diameter (mm)	length (cm)	approximate snare size (mm)	units per box
BX00711115	Exacto® cold snare	2.4	230	9	10

Distributed by:



5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639,4494

fax +1 440 / 639,4495

email global@usendoscopy.com

www.usendoscopy.com

¹ Horiuchi A, Hosoi K. "Prospective, randomized comparison of 2 methods of cold snare polypectomy for small colorectal polyps." *Gastrointestinal Endoscopy* (2015): 1-7.

² Din S, Ball A. "Cold snare polypectomy: Does snare type influence outcomes?" *Digestive Endoscopy* (2015): 1-6.
US Endoscopy is a subsidiary of STERIS Corporation.

Copyright© 2018 US Endoscopy or its affiliates. All rights reserved.

All marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

Follow us:    