

EUS balloons

Broadening our EUS product portfolio...

US Endoscopy is proud to offer a variety of latex and non-latex EUS balloon solutions.



The Ultramatrix™ EUS Balloons...

- are not made with natural rubber latex
- feature a polyisoprene material that provides a quality image
- are compatible with both Radial and Curvilinear echoendoscopes



Radial array



Curvilinear array



The Oracle® EUS non-latex balloon...

- is not made with natural rubber latex
- elastic band design allows for easy application and offers a secure grip to ensure the balloon stays on the endoscope
- is compatible with both Radial and Curvilinear echoendoscopes



Radial array



Curvilinear array

The Oracle® EUS latex balloons...

- provide a quality image
- are cost effective options
- are available for both Radial and Curvilinear echoendoscopes
- have been manufactured to contain latex



Radial array



Curvilinear array

Ultramatrix™ EUS balloons

product image	product number	description	sterile	units/box
	USBOL	Ultramatrix™ EUS balloon – Olympus linear*	no	20
	USBOR	Ultramatrix™ EUS balloon – Olympus / Pentax radial**	no	20

*Fits all Olympus linear GI EUS endoscopes.

** Fits all Olympus and Pentax radial GI EUS endoscopes.

These products have been cleared for use in endoscopic procedures that utilize a balloon in conjunction with an echoendoscope in the upper and lower gastrointestinal tract.

Oracle® EUS non-latex balloon

product image	product number	description	sterile	units/box
	BX00711891	Oracle® EUS non-latex balloon – Olympus curvilinear / radial array Fujifilm curvilinear*	yes	20

*The balloon is compatible with Olympus 140, 240, 160, 260 and Fujifilm curvilinear array echoendoscopes (models EG-530UT and EG-530UT2).

This product has been cleared for use in endoscopic procedures that utilize a balloon in conjunction with an echoendoscope in the upper and lower gastrointestinal tract.

Oracle® EUS latex balloons

product image	product number	package color	description	sterile	units/box
	BX00711892	navy blue	Oracle® EUS latex balloon – Olympus radial	yes	20
	BX00711893	baby blue	Oracle® EUS latex balloon – Olympus linear*	yes	20

*The balloons are compatible with Olympus 140, 240, 160, 260 and 180 curvilinear array echoendoscopes.

These products have been cleared for use in endoscopic procedures that utilize a balloon in conjunction with an echoendoscope in the upper and lower gastrointestinal tract.

For more information on our EUS balloon solutions, visit www.usendoscopy.com/eus.

Olympus® is a registered trademark of Olympus Optical, Inc.
 US Endoscopy is a subsidiary of STERIS Corporation.
 All marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.
 ©2018 US Endoscopy. All Rights Reserved.

Follow us:

5976 Heisley Road
 Mentor, OH 44060
 phone (1) 440 / 639-4494
 fax (1) 440 / 639-4495
 customer service (1) 800 / 769-8226
www.usendoscopy.com