

# Torrent™ irrigation pump & Torrent® irrigation system

**Torrent**  
irrigation pump



# A dependable total solution for intraprocedural irrigation...

Our full-line of irrigation solutions now includes the **new Torrent™ irrigation pump** - designed to provide consistent flow rates during endoscopic procedures when used in conjunction with the Torrent® irrigation tubing and Torrent® scope connector.

## The Torrent™ irrigation pump...

- offers a small, compact design for easy placement and transport on procedure carts and booms
- includes an adjustable flow rate knob allowing clinicians to increase/decrease flow based on type of procedure
- includes a replaceable pump head with a built-in safety feature - when opened, will not allow rollers to turn so the user can safely insert or remove tubing



Low-profile footswitch also available (EA00711565)

## Torrent® irrigation tubing & scope connectors

Proactively protect patients and staff against cross-contamination risks with the disposable 24-hour Torrent® irrigation tubing and single patient use Torrent® irrigation scope connector.

### *The Torrent® irrigation system...*

- features 24-hour use irrigation tubing – use for the day then throw it away
- offers a single patient use scope connector with a one-way luer-lock valve to help eliminate potential backflow



#### **Torrent® irrigation tubing**

Includes a 24-hour use red sticker to help identify date and time of setup



#### **Torrent® scope connector**

Single patient use, available for Olympus and Pentax endoscopes

The Food and Drug Administration (FDA) recently released an updated Guidance Document which recommends that scope connectors be discarded *“after every patient use.”*<sup>1</sup>

<sup>1</sup>Guidance For Industry and Food and Drug Administration Staff: “Mitigating the Risk of Cross-Contamination from Valves and Accessories Used for Irrigation through Flexible Gastrointestinal Endoscopes” November 2016.

## Product chart

product number	description	compatibility	units/box
BX00711569	Torrent™ irrigation pump The following accessories are included with each box: <ul style="list-style-type: none"> <li>• High profile footswitch (EA00711568)</li> <li>• Power cord</li> <li>• Water bottle holder (EA00711599)</li> <li>• Operating &amp; maintenance manual</li> </ul>	US Endoscopy Torrent® irrigation tubing (BX00711560/BX00711561) US Endoscopy Torrent® scope connectors (BX00711570/BX00711571)	1
BX00711560	Torrent® irrigation tubing (24-hour)	US Endoscopy Torrent™ irrigation pump Olympus OFF, OFF2, & AFU-100 pumps	25
BX00711561	Torrent® irrigation tubing (24-hour)	US Endoscopy Torrent™ irrigation pump Olympus OFF, OFF2, & AFU-100 pumps ERBE EIP2 pumps MEDIVATORS EndoGator & EndoStratus pumps	25
BX00711570	Torrent® scope connector (single patient use)	Olympus 140/160/180/190 series endoscopes	100
BX00711571	Torrent® scope connector (single patient use)	Pentax 70/80/90 series endoscopes	100
EA00711565	Low profile footswitch	US Endoscopy Torrent™ irrigation pump	1
EA00711568	High profile footswitch	US Endoscopy Torrent™ irrigation pump	1
EA00711599	Water bottle holder	US Endoscopy Torrent™ irrigation pump	1

**For more information on the Torrent™ irrigation pump and Torrent® irrigation system, contact US Endoscopy today.**

Olympus® is a registered trademark of Olympus Optical, Inc.  
ERBE USA Incorporated is a subsidiary of ERBE GmbH.

Medivators® is a registered trademark of Minrtach Corporation

PENTAX® is a registered trademark of Asahi Kogaku Kogyo Kabushiki Kaisha

US Endoscopy is a subsidiary of STERIS Corporation.

All marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy  
Copyright© 2018 US Endoscopy. All Rights Reserved.

5976 Heisley Road  
Mentor, OH 44060  
phone (1) 440 / 639-4494  
fax (1) 440 / 639-4495  
customer service (1) 800 / 769-8226

Follow us:    

[www.usendoscopy.com](http://www.usendoscopy.com)