








polypectomy snares



US endoscopy

listening... and delivering solutions®




From cold snaring to advanced polypectomy, US Endoscopy offers a product solution

	product number	description	needle gauge/ projection	sheath diameter (mm)	length (cm)	approximate snare size (cm)	units/ box
	BX00711018	Short Throw™ snare - mini	-	2.4	230	1.25 x 3.0	10
	BX00711019	Short Throw™ snare - standard oval	-	2.4	230	2.5 x 5.4	10
	BX00711101	Short Throw™ snare - standard oval enteroscopy	-	2.4	350	2.5 x 5.4	5
	BX00711084	iSnare® system - Lariat® snare	25 gauge x 5mm	2.5	230	oval: 30 x 50mm hex: 10 x 28mm diamond: 6 x 15mm	5
	BX00711085	iSnare® system - oval	23 gauge x 5mm	3.0	230	2.5 x 4.0	5
	BX00711086	iSnare® system - oval	25 gauge x 5mm	3.0	230	2.5 x 4.0	5
	BX00711088	iSnare® system - hexagonal	25 gauge x 5mm	3.0	230	2.5 x 4.0	5
	BX00711089	iSnare® system - hexagonal	23 gauge x 5mm	3.0	230	2.5 x 4.0	5
	BX00711113	Traxtion® hexagonal snare	-	2.4	230	2.5 x 4.0	10
	BX00711115	Exacto® cold snare	-	2.4	230	9.0mm	10

NEW!



on to fit a variety of clinical needs.

	product number	description	needle gauge/ projection	sheath diameter (mm)	length (cm)	approximate snare size (cm)	units/box
	BX00711117	Histolock® resection device	-	2.4	230	14mm x 27mm	10
	BX00711873	Histolock® resection kit - Histolock® resection device Carr-Locke injection needle Roth Net® retriever - polyp Mio® medical device organizer eTrap® polyp trap	-	-	-	-	3
	BX00711119	Lariat® snare	-	2.4	230	oval: 30 x 50mm hex: 10 x 28mm diamond: 6 x 15mm	10
	BX00711120	Rotator™ snare - standard	-	2.4	230	2.5 x 5.4	10
	BX00711121	Rotator™ snare - mini	-	2.4	230	1.25 x 3.0	10

All active cord connections are Olympus style.

For more information on the polypectomy snares, contact US Endoscopy today.



listening...and delivering solutions®

US Endoscopy is a subsidiary of STERIS Corporation.

©2018 US Endoscopy or its affiliates. All rights reserved.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

5976 Heisley Road

Mentor, OH 44060

phone (1) 440 / 639-4494

fax (1) 440 / 639-4495

customer service (1) 800 / 769-8226

www.usendoscopy.com

Follow us:    

760900 Rev. E