

CLEANING BRUSH REFERENCE CHART

	Cleaning Function	Product Code	Description	Channel Size	Brush Head	Bristle Diameter	Shaft Length	Qty/Box
	Valve	BX00711610	Valve/Control Head Cleaning Brush	-	Valve/Control Head	10mm	-	50
	Valve	BX00711611	Valve/Control Head Cleaning Brush	-	Valve/Control Head	10mm	-	100
	Channel & Valve	BX00711609	Double-Header™ combination channel cleaning brush & Valve/Control Head Cleaning Brush	2.0 - 7.0mm	Standard + Head Valve/Control	5-7mm taper	230cm	50
	Channel & Valve	BX00711629	Double-Header™ combination channel cleaning brush & Valve/Control Head Cleaning Brush	2.0 - 7.0mm	Standard + Head Valve/Control	5-7mm taper	230cm	100
	Channel & Valve	BX00711619	DuoSwift® combination squeegee brush + Valve/Control Head Cleaning Brush	2.8 - 7.0mm	Standard/Squeegee + Valve/Control Head	5-7mm (10mm squeegee)	230cm	50
	Channel	BX00711618	DuoSwift® Combination squeegee brush	2.8 - 7.0mm	Standard/Squeegee	5-7mm (10mm squeegee)	230cm	50
	Channel	BX00711601	Standard Channel Cleaning Brush	2.0 - 7.0mm	Standard	5-7mm taper	230cm	50
	Channel	BX00711607	Gemini Brush™ double-ended channel cleaning brush	2.0 - 7.0mm	Standard - Double	5-7mm taper	230cm	50
	Channel	BX00711627	Gemini Brush™ double-ended channel cleaning brush	2.0 - 7.0mm	Standard - Double	5-7mm taper	230cm	100
	Channel	BX00711612	Wooly Brush® channel cleaning brush	2.0 - 7.0mm	Standard - Length Increased	5-7mm taper	230cm	50
	Specialty Channel	BX00711603	Jumbo Channel Cleaning Brush	2.8 - 9.0mm	Standard - Bristle Diameter Increased	7-9mm taper	230cm	50
	Specialty Channel	BX00711606	Jumbo Channel Cleaning Brush	2.8 - 9.0mm	Standard - Bristle Diameter Increased	7-9mm taper	350cm	20
	Specialty Channel	BX00711616	Small Wonder® channel cleaning brush	1.0 - 2.25mm	Standard - Small Channel Specialty	2.25mm	160cm	15
	Specialty Channel	BX00711617	Small Wonder® channel cleaning brush	1.0 - 2.25mm	Standard - Small Channel Specialty	2.25mm	230cm	15