



BioGuard[®] air/water & suction valve kits air/water cleaning adapter



US Endoscopy is proud to expand our portfolio of disposable GI procedure product solutions. The new **BioGuard® air/water cleaning adapter - Olympus** joins the **BioGuard® air/water & suction valve kits** designed to help minimize cross-contamination concerns and support infection prevention practices.

The single-use BioGuard® air/water & suction valve kits offer an alternative option to manual cleaning and reprocessing and are available in 2, 3, and 4-piece configurations.

The BioGuard® air/water & suction valve kits...

- minimize the risk of cross-contamination between patients
- provide consistent performance – use a new kit for each procedure
- support adherence to industry guidelines which require traceability with a special endoscope¹



The new BioGuard® air/water cleaning adapter is designed specifically for Olympus GI endoscopes*, offering a single-use option that flushes the air/water channel during the bedside pre-cleaning process.

The BioGuard® air/water cleaning adapter...

- is disposable, offering an alternative option to manual reprocessing
- offers a unique appearance and warning tag to prevent use during patient procedures
- individually packaged to ensure use only during the bedside pre-cleaning process



For more information on the BioGuard® products, contact US Endoscopy today.

BioGuard® air/water suction valve kits			
product number	description	sterile	units/box
BX00711780	BioGuard® air/water & suction valves - 2 piece kit	yes	50
BX00711783	BioGuard® air/water & suction valves - 3 piece kit Kit includes: BioShield® biopsy valve (BX00711129)	yes	50
BX00711784	BioGuard® air/water & suction valves - 4 piece kit Kit includes: BioShield® biopsy valve (BX00711129) and Torrent® scope connector (BX00711570)	yes	50
<i>Each kit individually packaged and sealed in tray</i>			
BioGuard® air/water cleaning adapter - Olympus			
product number	description	sterile	units/box
BX00711788	BioGuard® air/water cleaning adapter - Olympus	no	50

*Compatible with Olympus GI endoscopes, excluding EUS echoendoscopes that utilize a balloon channel.

Distributed by:



5976 Heisley Road
Mentor, OH 44060 USA

phone: +1 440 / 639.4494
fax: +1 440 / 639.4495
email: global@usendoscopy.com
www.usendoscopy.com

¹ANSI/AAMI ST91: 2015 Flexible and semi-rigid endoscope processing in health care facilities. 20, 23, 25, 31, 32, 33.
Olympus® is a registered trademark of Olympus Corporation.
US Endoscopy is a subsidiary of STERIS Corporation.
©2019 US Endoscopy. All Rights Reserved.
All marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

