

# Centra™ disposable biopsy forceps



# An extensive portfolio of polypectomy devices to address a variety of clinical needs...

The cold Centra™ disposable biopsy forceps work with flexible endoscopes to obtain tissue samples for pathology analysis.

## The Centra™ disposable biopsy forceps feature a...

- variety of jaw designs - available with or without needle and serrated or plain cups
- distinct catheter color and position markers for visibility during insertion and withdrawal
- large jaw configuration opening over 90° for tissue sampling
- swinging jaw design to allow pivoting onto the targeted biopsy site

	product number	needle	jaw style	jaw OD	jaw opening (mm)	jaw length (mm)	length (cm)	minimum working channel (mm)	units/box
	BX00711245	no	plain cup	2.3	8	4	230	2.8	10
	BX00711247	yes	plain cup	2.3	8	4	230	2.8	10
	BX00711249	no	serrated	2.3	8	4	230	2.8	10
	BX00711250	yes	serrated	2.3	8	4	230	2.8	10

For more information on our specialty forceps available including hot, pediatric and enteroscopy, visit [www.usendoscopy.com](http://www.usendoscopy.com) today.



5976 Heisley Road  
Mentor OH 44060

phone 440 / 639,4494

fax 440 / 639,4495

customer service 800 / 769,8226

[www.usendoscopy.com](http://www.usendoscopy.com)

US Endoscopy is a subsidiary of STERIS Corporation.  
© 2019 US Endoscopy or its affiliates. All rights reserved.  
All marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

Follow us:    